



Rehab Car Seat **IPAI** **NXT**

... is flexible, comfortable and secure

5-point harness including quick lock adjustable backrest angle

grows with the child for a long time

a lot of innovative accessories



Rehab Car Seat **IPAI** **NXT**

...flexible, comfortable and secure!

The IPAI-NXT child seat is a for the disability field modified car seat. It is qualified for being installed in most of the customary vehicles. The car seat is designed to be used in the backseat or passenger seat of a car and requires a 3-point safety belt.

Additional to the car's seat belt, the IPAI-NXT comes with a 5-point-harness who gently protect your child. This can be adjusted in height together with the headrest to improve the belt route. By means of the continuous variable backrest angle, the seat can be adjusted to the car seat and to the child's needs. We installed soft cushion that improves the comfort enormously compared to a standard car seat.

Due to the headrest which is adjustable in many stages, the seat can be aligned to the size of the child. The IPAI-NXT provides an advanced lateral protection structure as he is equipped with fixed armrests in the hip area and energy absorption elements in the shoulder and chest area.

Due to the earated foam material, that is commonly used in the automotive industry, and an exclusive ventilation system that prevent from sweating, the childs optimal sitting comfort is ensured.





A popular feature is the innovated swivel which simplifies the transfer of the child even though it is gaining weight. The pivot point was moved from the center to a rear and outer position. This enables the seat to move away from the cars backrest during the rotation and simultaneously towards the doorway. During the development, great emphasis was placed on a compact design, so the total height increased just by 16 mm!

- 1 IPAI-NXT red, with swivel mechanism and adjustment knob
- 2 ISOFIX adapter, depth-adjustable to allow great adjustments of the seats and the backrest angles
- 3 continuous variable backrest angle (standard)
- 4 IPAI-NXT grey, with base board, angular adjustment 10°, swivel, hip guides and lateral supports



...multifarious, like your child's necessities

The IPAI-NXT is a child seat especially modified for children on rehabilitation and contains adaptable and individually positionable lateral supports to control the child's lateral movement. Optional accessories are available to ensure a comfortable and safe journey in the car:

- Angular adjustment (0 - 10° or 10 - 20°)
- Swivel
- Lateral supports
- ISOFIX adapter
- Security tray
- Abduction block
- Foot rest
- Extra cushion, 8-pieces
- Security pad for the buckle

By means of those accessories the individual needs of the child can be addressed, especially as the options are not mutually exclusive. Many modern vehicles are equipped with ISOFIX for which we provide suitable adapter units to install the seat perfectly in the car.

An important development is the 5-point-harness. As it can be tightened with only one move, it simplifies the daily life enormously.

We value your feedback and suggestions as we continuously develop the IPAI-NXT!



central quick lock for the 5-point-harness



gradual angular adjustment



headrest in many steps height-adjustable



foot rest adjustable in height and angle

Your supplier:

HERNIK GmbH
Bodenseestrasse 25
83059 Kolbermoor
☎ +49 (0)8031-2068044
☎ +49 (0)8031-2069215
✉ info@hernik.de
🌐 www.hernik.de