

invoice / delivery address

Reference:

☐ **Order**

☐ **Offer**

Rehab Car Seat **HERO** **NXT**

Exportprices

Release: 2024

**BASIC**

Details:

	HERO-NXT WITH STANDARD CUSHION	HERO-NXT WITHOUT STANDARD CUSHION
SEAT WIDTH (pelvis)	21	25 cm
SEAT WIDTH (knee)	28	28 cm
SEAT DEPTH	25	28 cm
BACK HEIGHT	21 (24) - 34	27 - 35 cm
CROWN HEIGHT	41 - 53	44 - 56 cm
DISTANCE LATERALS	12 - 21 cm	
total height	66 - 73 cm (without feet 60 cm)	
total width	43 cm	
total depth	63 cm	
capacity	18 kg	


common information (standard):

- crash tested according to ECE-R 129-00 (i.Size)
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- height adjustable back and head rest
- additional upholstery for the seat, back and head rest
- large seat inclination
- swivel 360° and adjustable stilt
- ISOFIX
- covers washable (30 °C)


E4

08025600  
UN Regulation  
No. 129/03

**Note:** ISOFIX adaptation in the car is absolutely necessary!!!



☐ **HERO-NXT  
red-black**  
(470000-E1)



☐ **HERO-NXT  
grey-black**  
(470000-E2)

### accessory for stabilization

abduction block

☐ (470400-E)




security tray

☐ (470600-E)



side supports made of  
aluminum fully adjustable

☐ (470800-E)



longer 5-point-harness  
(+10 cm)

☐ (470320-E)



## further accessories

protection mat  
(to protect original seats of the car)

☐ (470780-E)



head fixation  
(stabilizes the head if there's no control)

☐ (470630-E)



security pad for the buckle  
with snap fastener

☐ (470650-E)



belt to connect the chest straps

☐ (470640-E)



## notes